

LISTA DE INDICES DE REGRACCION / REFRACTIVE INDECES LIST

Medio de montaje / <i>Embedding medium</i>	Indice de refraccion / <i>Refractive index</i>
Agarosa / <i>Agarose</i>	1.34
Fluoromount-G	1.40
Medio de Cultivo / <i>Growth medium</i> (Experimentos in vivo / <i>In vivo</i> experiments)	1.34
Glicerol / <i>Glycerol</i> (90%)*	1.46
Glicerol / <i>Glycerol</i> (75%)	1.44
Glicerol / <i>Glycerol</i> (50%)**	1.40
Mowiol	1.41 – 1.49
ProLong Diamond	1.47
ProLong Glass	1.52
Vectashield	1.45
Vectashield HardSet (dry)	1.46
Vectashield Vibrance	1.47
Agua / <i>Water</i>	1.33
Indice de refraccion de la lente (medio de inmersión) <i>Lens refractive index (immersion medium)</i>	Indice de refraccion / <i>Refractive index</i>
Aire / <i>Air</i>	1.00
Glicerol / <i>Glycerol</i> (90%)*	1.46
Glicerol / <i>Glycerol</i> (75%)	1.44
Glicerol / <i>Glycerol</i> (50%)**	1.40
Agua / <i>Water</i>	1.333
Zeiss Immersol 518F	1.518
Zeiss Immersol W	1.334
Olympus Immersion Oil F	1.518
Olympus Silicone Immersion Oil	1.406

* Se usa solo para fluorescencia / Used for fluorescence images only

** Se usa solo para imágenes DIC / Used for DIC images only

29/09/2020